

Multi-function Home Care Bed

(Back-rest and Knee-rest simultaneously)

(Comply to Safety Standard)



Remote controller

EN-1M

We pursued safety and functionality!

Movable Fences Siderail

Easy attached & detached while using




Remote controller

Allowing control of back, knee function




Movable Armrest Siderail

Armrest with non-slip cover could spin out and grab firmly to support patient while standing up.






Full-length Siderails

It made of ABS plastic and provides the optimum protection against a fall.




Features



- The EN-1M Bed's frame is made of steel with powder coating
- The bed surface is made of steel with mesh design; Head & Foot board are made of ABS plastic with plinth design, providing better stability.
- It is easy to set up the bed height; there are 20, 24, 28, 32, 36, 40cm for choosing.
- This function manipulates the back-rest and knee-rest simultaneously to the position of a chair; it provides a comfortable position for user.
- Total weight is 50kg, it is easy to move and delivery.

Packaging Dimensions


1. Head and foot board : 102 x 69 x 11 cm
2. Beams: 206 x 17 x 9 cm
3. Bottom frame: 93 x 88 x 16 cm
4. Folding bed frame: 93 x 88 x 7 cm




Specifications




Back Angle 0° ; Knee Angle 0°



Back Angle 13° ; Knee Angle 15°



Back Angle 45° ; Knee Angle 12°



Back Angle 75° ; Knee Angle 0°

