



# ES-02FDS MS-02FS MF-02FS Electric & Manual Hospital Bed




Collapsible Cranks, for manually operating back, knee & Hi-Lo functions.



Two gas springs combine with Hi-Lo rod can extend life of crank.



Allowing control of back, knee, Hi-Lo & auto-contour functions.



Head/Foot Board Lock, for securing head/foot board.



5" Dust-Proof Central Braking Module, operating by pedal lever located at foot-end of the bed for one step brake.

**Movable Third Siderail**  
Optional feature to add on bed for upgraded protection.



**Attendant Control Panel**  
This provides caregivers easy control of the position of the patient.

**Lock-Out Function**

Allowing lock-out all controls of bed to prevent unexpected touch.



**Auto-Regression**

Release all positions and height to be flat position.