



ES-96HD

Electric Hospital Bed




X-Ray Cassette



Backrest X-Ray Module, allowing examination without patient mobilization.



Siderail Controls, inner control allowing patient operates back, knee, Hi-Lo & auto-contour functions; and outer control provides reading of backrest angle & prevent infection.



Siderail Control





5" Dust-Proof Central Braking & Steering Module, operating by pedal lever located at foot-end of the bed for one step brake.



CPR and Quick Release

Backrest CPR Lever, manually releasing backrest for CPR operation.



ACP includes Lock-Out, Auto-Regression controls.