JA-107

Emergency Cart

JOSON-CARE ENTERPRISE CO., LTD

Http://www.hospitalbed-josoncare.com/ E-mail:overseas_ma@joson-care.com.tw Factory:

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Township, Taoyuan County 33348, Taiwan (R.O.C.)

Tel:886-3-3290925 Fax:886-3-3290921

No.280, Chengdu Rd., Xitun Dist., Taichung City 40751, Taiwan(R.O.C.)

Tel:886-4-23160743 Fax:886-4-23160740





Two options on frame coating colors:



Blue



Optional Accessories



Features:

Orange

- 1. The frame of cart is made by steel with epoxy powder coating.
- 2. The upper and lower covers of frame are made by ABS for easy clean and attraction.

Lock Piece with Disposable Lock

- 3. Drawers have trigger design to prevent sliding out while moving.
- 4. Security lock function is optional manipulated by either key or disposable lock.
- 5. Drawers include 3" x 2 sets, 6" x 2 sets and 9" x 1 set.
- 6. There are two clapboards of each vertical and horizontal in both 3" and 6" drawers.
- 7. Diversity accessories with easy set up and flexible configuration design.





Holder







Glove Box

Glove Box Holder