



Safety Quality Sustainability



Business
Concept →

Mission

Serving all mankind through providing professional medical products

Vision

Improving healthcare environment, quality and efficiency

Core Value

- Safety
- Quality
- Sustainability

History →

2000

Established

Joson-Care Enterprise Co., Ltd.

2007



Compliance

IEC 60601-1
IEC 60601-1-2

2005



Certificate Approved

ISO 13485, ISO 9001, CE & FDA

Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.



FDA



TFDA



MDR



IEC60601-1




IEC60601-2



IEC60601-2-38



IEC60601-2-52



So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.

2008



Compliance
60601-2-38

2014




Compliance
60601-2-52

2023



MDR

Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.



ICU Bed

5



- 19 series
- ES-19HD ICU Bed 7
 - ES-19DF ICU Bed 8



- 12 series
- ES-12DW ICU Bed with Weight Scale 13
 - ES-12HD Electric Hospital Bed 15
 - ES-12DF Electric Hospital Bed 16
 - ES-11DF Electric Hospital Bed 17



- 09.96 series
- ES-09DW Electric Hospital Bed with Weight Scale 21
 - ES-96HD Electric Hospital Bed 21

Hospital Bed

23



- ES-05HD Electric Hospital Bed 25
- ES-05HDS Electric Hospital Bed 26
- MS-05HS Manual Hospital Bed 27
- ES-99HDS Electric Hospital Bed 27



- 02.16 series
- ES-02JH Electric Nursing Bed 31
 - ES-16N series Electric Nursing Bed 32
 - MS-16NS Manual Hospital Bed 33

Emergency Stretcher

35



- Stretcher series
- JE-350 & JE-300 Hydraulic Stretcher 37
 - JE-200 & JE-150 Hydraulic Stretcher 39
 - JE-130 3 Manual Cranks Emergency Stretcher 42
 - JE-120 2 Manual Cranks Emergency Stretcher 42
 - JE-100 Simple Emergency Stretcher 42

Pediatric Bed

43



Pediatric Bed series

● JC-012 Electric Pediatric Bed	45
● JC-011 Manual Pediatric Bed	46
● JC-031 Infant Crib	47
● JC-030 Crib	47

Medical & Ward Equipment

Mattress	67
JT-100 Orthopedic Traction Frame	69
IV Stand	70
JA-107 Emergency Cart & UD Medical Cart	71
Medical Cart	73
Cabinet & Outpatient Service Equipment	74
Overbed Table	75
Caregiver Recliner & Doctor Chair	76

Home Care Bed

49



18 series

● ES-18UDS Ultra Low Bed	51
● ES-17UDS Ultra Low Bed	52



EN series

● EN-3M Foldable Electric Home Care Bed	57
● EN-2M Assembly Electric Home Care Bed	58



08 series

● ES-08FDS Electric Hospital Bed	63
● ES-08FSS Electric Hospital Bed(Metal frame)	64
● MS-08FS 3 Manual Cranks Hospital Bed	65
● MF-08FS 2 Manual Cranks Hospital Bed	65

19 series




ES-19HD / ES-19DF



Head Elevation
0 ~ 75°



Knee Elevation
0 ~ 30°



Bed Elevation



Auto-chair



Trend. $\pm 18^\circ$



Rev. Trend. $\pm 18^\circ$



Trend.
Angle Indicator




Head
Angle Indicator




Safety Standard



230KG
Work Load



4PCS
Actuator



Backup Battery



IP66
Actuator/Control
Box





Shear Stress
Bi-directional Release

ES-19HD

ICU Bed

Standards: IEC 60601-1, 60601-1-2, 60601-2-52



Shear Stress
Bi-directional Release



Bed Elevation 40cm



Backrest X-Ray
Module



Trend. / Rev. Trend. $\pm 18^\circ$



Head Angle Indicator /
Outer Control Panel





Night Light


ES-19DF

ICU Bed

Standards: IEC 60601-1, 60601-1-2, 60601-2-52



Chair Position



Bed Height Indicator



Brake Light Indicator



Storage Space



Bed-end Extension 20cm




Mattress Stopper


19 series

Shear Stress Bi-directional Release

The hip and back boards decline to minimize shear stress and ease pressure on the back and abdomen.



Bed Elevation 40cm



LO: 38 cm

● Green: In low pos.

HI: 78cm


● Orange: Above 45cm

Vertical Lifting

The bed's frame is fixed while elevation is being operate.



Non-moving Headboard



Headboard fixed during Trend, and elevation adjust.



Specification	ES-19HD	ES-19DF
Overall bed length	L : 2220mm	L : 2220mm
Overall bed width	W : 958mm	W : 958mm
(Siderails down/up)	990mm	990mm
Overall mattress base length	L : 1925mm	L : 1925mm
Overall mattress base width	W : 885mm	W : 885mm
Overall bed height	(5") H : 385~785mm (6") H : 410~810mm	(5") H : 385~785mm (6") H : 410~810mm
Head Elevation	0~75°	0~75°
Knee Elevation	0~30°	0~30°
Trend./Rev. Trend.	±18°	±18°
Work Load	230kg	230kg
Castor Diameter	5" or 6"	5" or 6"
Actuator	4	4
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66
Safety Standard	IEC 60601-2-52	IEC 60601-2-52

Applicable mattress size (mm)	L1900xW890xH100	L1900xW890xH100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Features	ES-19HD	ES-19DF
Backrest X-Ray Module	☒	☐
Caregiver Siderail Control	☒	☒
Patient Siderail Control	☒	☒
LCD Screen Panel Remote Controller	☒	☐
Backrest, Trend./Rev. Trend. Angle Indicator	☒	☒
CPR Quick Release	☒	☐
Roller Bumper	☒	☒
Lock-Out Function	☒	☒
Auto Regression	☒	☒
Examination Position	☒	☒
Anti-Shock Position	☒	☒
Chair Position	☒	☒
Trend./Rev. Trend.	☒	☒
Bed-end Extension	☒	☒
Dust-Proof Central Braking/Steering Module	☒	☒
IV Pole Hole	☒	☒
Drainage Bag Hanger	☒	☒
Mattress Stopper	☒	☒
Brake/Height Caution Indicator	☒	☒
Night Light	☒	☐
Third Siderail	☐	☐
Restraint Belt Hole	☐	☐

☒ Standard ☐ Option


12 series



ES-12HD / ES-12DW




ES-12DF / ES-11DF



ES-12DW

ICU Bed
with Weight Scale

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



Weight Monitoring Module



Bed Height Indicator



Trend. / Rev. Trend. $\pm 12^\circ$



Chair Position



Mattress Stopper

11 / 12 series


Safety Standard

11、12 series comply to latest safety standard IEC 60601-2-52, which ensures that all gaps between head/foot boards and siderails are smaller than 60mm



ES-12DW

Weight Scale & ACP Functions



Fall Prevention Alarm





Weight Sensors

Sensors placed around the bed to detect when a patient is attempting to get out of bed

ES-12HD

Electric Hospital Bed

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



5-function ACP including
Lock-Out and Auto-Regression control



X-Ray Cassette



Head Angle Indicator /
Trend. Angle Indicator





Bed-end Extension 20cm

ES-12DF

Electric Hospital Bed

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



Chair Position



Bed Height Indicator



Safety Standard
IEC 60601-2-52





Bed-end Extension 20cm

ES-11DF

Electric Hospital Bed

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



Trend. / Rev. Trend. $\pm 12^\circ$
(Option)



CPR Quick Release
(Option)



Drainage Bag Hanger



Bed-end Extension 10cm



Mattress Stopper



Specification	ES-12DW	ES-12HD	ES-12DF	ES-11DF
Overall bed length	L : 2220mm	L : 2220mm	L : 2220mm	L : 2120mm
Overall bed width (Siderails down/up)	W : 950mm 1000mm	W : 950mm 1000mm	W : 950mm 1000mm	W : 950mm 1000mm
Overall mattress base length	L : 2020mm	L : 2020mm	L : 2020mm	L : 1920mm
Overall mattress base width	W : 900mm	W : 900mm	W : 900mm	W : 900mm
Overall bed height	H : 455~755mm	H : 455~755mm	H : 455~755mm	H : 455~755mm
Head Elevation	0~75°	0~75°	0~75°	0~75°
Knee Elevation	0~30°	0~30°	0~30°	0~30°
Trend./Rev. Trend.	±12°	±12°	Optional	Optional
Work Load	230kg	230kg	230kg	230kg
Castor Diameter	5"	5"	5"	5"
Actuator	4	4	3/4	3
Volts	AC 100-240V	AC 100-240V	AC 100-240V	AC 100-240V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah	DC 24V 2.0Ah	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66	IP66	IP66
Safety Standard	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52

Applicable mattress size (mm)	L2000xW880xH120	L2000xW880xH120	L2000xW880xH120	L1900xW880xH120
-------------------------------	-----------------	-----------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

Features	ES-12DW	ES-12HD	ES-12DF	ES-11DF
X-Ray Cassette	✓	✓	□	□
Outer/Inner Control Panel	✓	✓	✓	✓
Head Angle Indicator	✓	✓	✓	✓
Trend./Rev. Trend. Angle Indicator	✓	✓	✓	✓
CPR Quick Release	✓	✓	□	□
Roller Bumper	✓	✓	✓	✓
Attendant Control Panel	✓	✓	✓	✓
Remote Controller	□	□	□	□
Weight Scale	✓			
Lock-Out Function	✓	✓	✓	✓
Auto Regression	✓	✓	✓	✓
Backup Manual Control Module			□	□
Bed-end Extension	✓	✓	✓	✓
Chair Position	✓	✓	✓	
Trend./Rev. Trend.	✓	✓	□	□
Dust-Proof Central Braking/Steering Module	✓	✓	✓	✓
IV Pole Hole	✓	✓	✓	✓
Drainage Bags Hanger	✓	✓	✓	✓
Mattress Stopper	✓	✓	✓	✓
Height Caution Indicator	✓	✓	✓	✓
Night Light	□	□	□	□
LCD Screen Panel Remote Controller	□	□	□	□

✓ Standard □ Option

09.96 series



ES-09DW




Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 22°



Bed Elevation



Auto-chair



Trend. $\pm 12^\circ$



Rev. Trend. $\pm 12^\circ$



Head Elevation
indicator




Safety Standard




225KG

Work Load



4PCS

Actuator



Backup Battery



IP66

Actuator/Control
Box

ES-96HD




Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°



Bed Elevation



Auto-chair



Trend. $\pm 12^\circ$



Rev. Trend. $\pm 12^\circ$



Trend.
Angle Indicator



Head
Angle Indicator




Safety Standard




225KG

Work Load



4PCS

Actuator



Backup Battery





IP66

Actuator/Control
Box

ES-09DW

Electric Hospital Bed
with Weight Scale




Standards: IEC 60601-1-2



ES-96HD

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Attendant
Control Panel (ACP)



Drop Down Siderail



Mattress Stopper



Specification	ES-09DW	ES-96HD
Overall bed length	L : 2250mm	L : 2250mm
Overall bed width	W : 940mm	W : 940mm
(Siderails down/up)	1060mm	1060mm
Overall mattress base length	L : 2020mm	L : 2020mm
Overall mattress base width	W : 925mm	W : 925mm
Overall bed height	H : 455~755mm	H : 455~755mm
Head Elevation	0~85°	0~85°
Knee Elevation	0~22°	0~32°
Trend. / Rev. Trend.	±12°	±12°
Work Load	225kg	225kg
Castor Diameter	5"	5"
Actuator	4	4
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66
Safety Standard	IEC 60601-2-38	IEC 60601-2-38
Applicable mattress size (mm)	L1980×W890×H120	L1980×W890×H100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Features	ES-09DW	ES-96HD
X-Ray Cassette	☒	
Outer/Inner Control Panel	☒	☒
Backrest, Trend./Rev. Trend. Angle Indicator	☒	☒
CPR Quick Release	☒	☐
Bumper	☒	☒
Attendant Control Panel	☒	☒
Electric CPR	☒	☒
Remote Controller	☐	☐
Weight Scale	☒	
Lock-Out Function	☒	☒
Auto Regression	☒	☒
Lockable Castor		☒
Dust-Proof Central Braking/Steering Module	☒	☐
Dust-Proof Lockable Caster		☐
IV Pole Hole	☒	☒
Drainage Bag Hanger	☒	☒
Mattress Stopper	☒	☒

☒ Standard ☐ Option

05.99 series




ES-05HD / ES-05HDS




Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°



Bed Elevation




Auto-chair



Safety Standard



225KG
Work Load



3PCS
Actuator




Backup Battery



IP66

Actuator/Control
Box


MS-05HS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°



Bed Elevation




225KG

Work Load



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°




Bed Elevation



225KG

Work Load



3PCS

Actuator



Backup Battery



Auto-chair



Safety Standard



IP66


Actuator/Control
Box

ES-99HDS

ES-05HD

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Attendant Control Panel (ACP)



Swivel Siderail



Safety Design for Handling Siderail




Backrest Angle Indicator (Option)

ES-05HDS


Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Siderails up

Bed width 100cm fit to small ward.



100cm



Detachable
Head/Foot Board



IV Pole Hole





Drainage Bag Hanger



Mattress Stopper

MS-05HS



Manual Hospital Bed




ES-99HDS

Electric Hospital Bed

Standard: IEC 60601-1-2, IEC 60601-2-38



Manual Crank Screw drive



Mattress Stopper



Swivel Siderail



IV Pole Hole



Restraint Belt Hole



Specification	ES-05HD	ES-05HDS	MS-05HS	ES-99HDS
Overall bed length	L : 2170mm	L : 2115mm	L : 2170mm	L : 2115mm
Overall bed width (Siderails down/up)	W : 940mm 1090mm	W : 910mm 1000mm	W : 940mm 1090mm	W : 930mm 1090mm
Overall mattress base length	L : 1940mm	L : 1940mm	L : 1940mm	L : 1940mm
Overall mattress base width	W : 925mm	W : 845mm	W : 925mm	W : 925mm
Overall bed height (Central Braking/Normal Castor)	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm
Head Elevation	0~85°	0~85°	0~85°	0~85°
Knee Elevation	0~32°	0~32°	0~32°	0~32°
Trend./Rev. Trend.	None	None	None	None
Work Load	225kg	225kg	225kg	225kg
Caster Diameter	5"	5"	5"	5"
Actuator	3	3	Manual Crank	3
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz	None	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah	None	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66	None	IP66
Safety Standard	IEC 60601-2-38	IEC 60601-2-38	None	IEC 60601-2-38

Applicable mattress size (mm)	L1900xW890xH100	L1930xW840xH100	L1900xW890xH100	L1900xW890xH100
-------------------------------	-----------------	-----------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree




Features	ES-05HD	ES-05HDS	MS-05HS	ES-99HDS
Outer/Inner Siderail Control	✓	✓		✓
Backrest, Trend./Rev. Trend. Angle Indicator	□	□	□	□
CPR Quick Release	□	□		□
Bumper	✓		✓	
Attendant Control Panel	□	□		
Electric CPR	✓			
Remote Controller	□	□		□
Lock-Out Function	□	□		
Auto Regression	✓			
Pneumatic Cylinder assisted Hi-Lo Mechanism			✓	
Manual Trend./Rev. Trend.	□	□	□	□
Dust-Proof Central Braking/Steering Module	□	□	□	□
Dust-Proof Lockable Castor	□	□	□	□
Lockable Castor	✓	✓	✓	✓
IV Pole Hole	✓	✓	✓	✓
Drainage Bag Hanger	✓	✓	✓	✓
Mattress Stopper	✓	✓	✓	✓

✓ Standard □ Option



02.16 series






ES-02JH

								
Head Elevation 0 ~ 85°	Knee Elevation 0 ~ 32°	Bed Elevation	Auto-chair	Safety Standard	170KG Work Load	3PCS Actuator	Backup Battery (Option)	IP54 Actuator/Control Box

ES-16N

							
Head Elevation 0 ~ 85°	Knee Elevation 0 ~ 40°	Bed Elevation	Safety Standard	225KG Work Load	3PCS Actuator	Backup Battery (Option)	IP54 Actuator/Control Box


MS-16NS

				
Head Elevation 0 ~ 85°	Knee Elevation 0 ~ 40°	Bed Elevation	Safety Standard	225KG Work Load

ES-02JH

Electric Nursing Bed

Standards: IEC 60601-1-2, IEC 60601-2-52



Rotary Armrest Siderail (Option)



Locking mechanism for fixed angles.



Detachable Head/Foot Board



Siderail Holder



Remote Controller



IV Pole Hole




Mattress Stopper

ES-16N series

Electric Nursing Bed

Standards: IEC 60601-1-2, IEC 60601-2-52

3/4 Length Siderail




3/4 Length Siderail + P-Type Siderail




Attendant Control Panel
(Option)



Remote Controller



Detachable
Head/Foot Board




Mattress Stopper


MS-16NS

Manual Hospital Bed


Standards: IEC 60601-1-2, IEC 60601-2-52



3/4 Length Siderail + P-Type Siderail




Manual Screw drive



IV Pole Hole



Detachable
Head/Foot Board



Mattress Stopper



Specification	ES-02JH	ES-16N	MS-16NS
Overall bed length	L : 2130mm	L : 2115mm	L : 2115mm
Overall bed width	W : 1000mm	W : 1000mm	W : 1000mm
Overall mattress base length	L : 1910mm	L : 1910mm	L : 1910mm
Overall mattress base width	W : 845mm	W : 845mm	W : 845mm
Overall bed height	H : 320~630mm	H : 350~660mm	H : 350~660mm
Head Elevation	0~75°	0~85°	0~85°
Knee Elevation	0~30°	0~40°	0~40°
Work Load	170kg	225kg	225kg
Castor Diameter	4"	5"	5"
Actuator	3	3	Manual Crank
Volts	AC 120V	AC 120V	None
Frequency	50/60Hz	50/60Hz	None
Backup Battery	Optional	Optional	None
Actuator/Control Box	IP54	IP54	None
Safety Standard	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52

Features	ES-02JH	ES-16N	MS-16NS
Remote Controller	☒	☐	
Siderail	☐	☒	☒
Bumper	☒	☒	☒
Attendant Control Panel	☐	☐	
Lock-Out Function		☐	
Auto Regression		☐	
Pneumatic Cylinder assisted Hi-Lo Mechanism			☒
Dust-Proof Central Braking /Steering Module	☐	☐	☐
Dust-Proof Lockable Castor	☒	☐	☐
Lockable Castor		☒	☒
IV Pole Hole	☒	☒	☒
Drainage Bag Hanger	☒	☒	☒
Mattress Stopper	☒	☒	☒

☒ Standard ☐ Option

Applicable mattress size (mm)	L1930xW840xH100	L1930xW840xH100	L1930xW840xH100
-------------------------------	-----------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree (ES-02JH ±3°)

Stretcher series



JE-350 & JE-300



Head Elevation
0 ~ 90°



Bed Elevation



Trend. $\pm 16^\circ$




Rev. Trend. $\pm 16^\circ$



225KG
Work Load

JE-200 & JE-150




Head Elevation
0 ~ 90°



Knee Elevation
0 ~ 30°



Bed Elevation



Trend. $\pm 18^\circ$



Rev. Trend. $\pm 18^\circ$






225KG
Work Load

JE-350


Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



unit: mm



X-Ray Cassette



Whole Body X-Ray
Cassette System



Functional Pedal






Monitor Tray

JE-300


Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



unit: mm



Back Gas Spring



Siderail Lock Handgrip



Oxygen Cylinder



3-section IV Pole






Foot Gas Spring
(Option)

JE-200

Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



unit: mm



X-Ray Cassette
(Option)



Retractable Push Handgrip





Oxygen Cylinder



Paper Roll Holder

JE-150

Hydraulic Stretcher



Lifting Rod



Knee Elevation Crank



Collapsible Siderail



Siderail Lock Handgrip



Dust-Proof Central Braking
Steering Module



Specification	JE-350	JE-300	JE-200	JE-150	Features	JE-350	JE-300	JE-200	JE-150
Overall bed length	L : 2120mm	L : 2120mm	L : 2100mm	L : 2100mm	X-Ray	✓		□	
Overall bed width	W : 790mm	W : 790mm	W : 825mm	W : 825mm	Full-Length X-Ray Capacity	✓			
(Siderails down/up)	W : 830mm	W : 830mm	W : 840mm	W : 840mm	Outward Siderail	✓	✓		
Overall mattress base length	L : 1850mm	L : 1850mm	L : 1900mm	L : 1900mm	Collapsible Siderail			✓	✓
Overall mattress base width	W : 635mm	W : 635mm	W : 740mm	W : 740mm	Knee Elevation Handgrip	✓	✓	✓	✓
Overall bed height	H : 590~890mm	H : 590~890mm	H : 585~935mm	H : 585~935mm	Roller Bumper	✓	✓		
Head Elevation	0~90°	0~90°	0~90°	0~90°	Bumper			✓	✓
Knee Elevation	None	None	0~22°	0~22°	Dust-Proof Central Braking/Steering Module	✓	✓	✓	✓
Trend./Rev. Trend.	16°	16°	18°	18°	The Fifth Wheel Module	✓	✓	✓	
Work Load	225kg	225kg	225kg	225kg	Oxygen Bottle Holder	✓	✓	□	□
Frame Material	Power Coated Mild Steel	Power Coated Mild Steel	Power Coated Mild Steel	Power Coated Mild Steel	Storage Tray	✓	✓	✓	✓
Castor Diameter	8"	8"	8"	8"	IV Pole	✓	✓	✓	✓
Transmission	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held; Crank Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held; Crank Bed Hi-Lo: Hydraulic Cylinder	Drainage Bag Hanger	✓	✓	□	□
					Restraint Belt Hole	□	□	□	□
Applicable mattress size(mm)	L1820xW620xH75	L1820xW620xH75	L1850xW720xH100	L1850xW720xH100	✓ Standard □ Option				

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

JE-130

3 Manual Cranks
Emergency Stretcher



JE-120

2 Manual Cranks
Emergency Stretcher



JE-100

Emergency Stretcher



Specification	JE-130	JE-120	JE-100
Overall bed length	L : 1900mm	L : 1900mm	L : 1900mm
Overall bed width	W : 750mm	W : 750mm	W : 750mm
Overall mattress base length	L : 1810mm	L : 1810mm	L : 1810mm
Overall mattress base width	W : 550mm	W : 550mm	W : 550mm
Overall bed height	H : 550~750mm	H : 650mm	H : 650mm
Head Elevation	0~80°	0~80°	0~80°
Knee Elevation	0~35°	0~35°	None
Work Load	150kg	150kg	150kg
Frame Material	SUS304	SUS304	SUS304
Castor Diameter	6"	6"	6"
Transmission	Crank	Crank	Hand-Hold
Mattress	2" High density PU foam	2" High density PU foam	2" High density PU foam

Applicable mattress size (mm)	L1800xW550xH50	L1800xW550xH50	L1800xW550xH50
-------------------------------	----------------	----------------	----------------

±20mm tolerance in all dimensions

±30mm in siderail width

+5° -0° tolerance in angle degree

Features	JE-130	JE-120	JE-100
Collapsible Siderail	✓	✓	✓
Lockable Castor	✓	✓	✓
Oxygen Bottle Holder	✓	✓	✓
Storage Tray	✓	✓	✓
IV Pole	✓	✓	✓

✓ Standard

Pediatric Bed series



JC-012



Head Elevation
0 ~ 80°



Knee Elevation
0 ~ 40°




Bed Elevation




170KG

Work Load



3PCS

Actuator



Backup Battery



IP66

Actuator/Control
Box


JC-011



Head Elevation
0 ~ 80°



Knee Elevation
0 ~ 40°




90KG

Work Load



Head Elevation
0 ~ 80°




90KG

Work Load

JC-010

JC-012

Electric Pediatric Bed



Safety Siderail Handle

Attendant Control Panel (ACP)



Attendant Control Panel (ACP)



Safety Siderail Handle




Remote Controller



Dust-Proof Lockable Castor

JC-011

Manual Pediatric Bed



Manual 2 Cranks




Solid Structure



Safety Siderail Handle

JC-010

Manual Pediatric Bed



JC-031 Infant Crib

Size : L700 x W440 x H910mm
 ± 20 mm tolerance in all dimensions
 Frame : ABS Plastic, Acrylic Basin
 Basin Angle: 0~20°



JC-030 Infant Crib

Size : L680 x W380 x H760mm
 ± 20 mm tolerance in all dimensions
 Frame : SUS304, Acrylic Baby Basin



Specification	JC-012	JC-011	JC-010
Overall bed length	L : 2050mm	L : 2000mm	L : 1600mm
Overall bed width	W : 1000mm	W : 816mm	W : 816mm
Overall mattress base length	L : 1870mm	L : 1944mm	L : 1545mm
Overall mattress base width	W : 845mm	W : 756mm	W : 758mm
Overall bed height	H : 650~850mm	H : 763mm	H : 764mm
Head Elevation	0~80°	0~80°	0~80°
Knee Elevation	0~40°	0~40°	None
Work Load	170kg	90kg	90kg
Frame Material	SUS304+Powder Coated	SUS304+Powder Coated	SUS304+Powder Coated
Castor Diameter	5"	4"	4"
Actuator	3	Manual (2 Cranks)	Manual (1 Crank)
Volts & Frequency	AC 100-240V ; 50/60Hz	None	None
Backup Battery	DC 24V 2.0Ah	None	None
Actuator/Control Box	IP66	None	None

Applicable mattress size (mm)	L1850xW830xH30	L1920xW720xH30	L1440xW780xH76
-------------------------------	----------------	----------------	----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Features	JC-012	JC-011	JC-010
Remote Controller	☒		
Bumper	☒		
Attendant Control Panel	☒		
Electric CPR	☒		
Lock-Out Function	☒		
Auto Regression	☒		
Dust-Proof Central Braking/Steering Module	☒		
Dust-Proof Lockable Castor	☐		
Lockable Castor		☒	☒
IV Pole Hole	☒	☒	☒

☒ Standard ☐ Option

18 series



ES-18UDS / ES-17UDS



Head Elevation
0 ~ 70°



Knee Elevation
0 ~ 30°




Bed Elevation



Auto-chair



170KG
Work Load



3PCS
Actuator



Backup
Battery Box
(Option)





IP54
Actuator/Control
Box




Shear Stress
Bi-directional Release

ES-18UDS

Ultra Low Bed



LO: 23.5cm



Shear Stress
Bi-directional Release



Dust-Proof Central Braking
Steering Module





Detachable
Head/Foot Board



Mattress Stopper

ES-17UDS

Ultra Low Bed



Shear Stress
Bi-directional Release



Detachable
Head/Foot Board



Dust-Proof Central Braking
Steering Module




3/4 Length Siderail, Rotary Armrest Siderail, Plug-in Siderail

18 / 17 series

Ultra Low Design

Ultra low at **23.5 cm**



ES-18UDS

23.5cm



ES-17UDS

Shear Stress Bi-directional Release

The hip and back boards decline to minimize shear stress and ease pressure on the back and abdomen.



More space between back board and hip board.





Specification	ES-18UDS	ES-17UDS
Overall bed length	L : 2080mm	L : 2100mm
Overall bed width	W : 1050mm	W : 950mm
Overall mattress base length	L : 1915mm	L : 1915mm
Overall mattress base width	W : 910mm	W : 830mm
Overall bed height	4"Lockable Caster 235~535mm	4"Lockable Caster 235~535mm
Head Elevation	0~70°	0~70°
Knee Elevation	0~30°	0~30°
Work Load	170kg	170kg
Castor Diameter	4"	4"
Actuator	3	3
Volts & Frequency	AC110V ; 50/60Hz	AC110V ; 50/60Hz
Actuator/Control Box	IP54	IP54

Applicable mattress size (mm)	L1900×W890×H100	L1830×W830×H100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Features	ES-18UDS	ES-17UDS
Remote Controller	☒	☒
Rotary Siderail	☐	☐
3/4 Length Siderail	☐	☐
Plug-in Siderail	☐	☐
Detachable Head/Foot Board	☒	☒
Dust-Proof Central Braking / Steering Module	☒	☒
Dust-Proof Lockable Castor	☐	☐
Lockable Castor	☐	☐
IV Pole Hole	☒	☒
Mattress Stopper	☒	
Backup Battery (option) 9V Alkaline Battery x3	☐	☐

☒ Standard ☐ Option


EN series



EN-3M






EN-2M



EN-3M

Foldable Electric Home Care Bed



Trend./Rev. Trend. $\pm 11^\circ$



Chair Position



Wooden Head/Foot Board





3/4 Length Siderail, Rotary Armrest Siderail, Plug-in Siderail

EN-2M

Assembly Electric Home Care Bed

5 Cartons Pack



Trend./Rev. Trend. $\pm 11^\circ$



Chair Position



Wooden Head/Foot Board



Dust-Proof
Lockable Castor



Mattress Stopper

3 Siderails




Rotary Armrest Siderail

Siderail installs on bedside, caregivers could adjust angles upon patient's need.



Plug-in Siderail

Siderail installs on bedside or foot-end.



3/4 Length Siderail

Metal tenon used

Siderails for below models

ES-02JH | EN-3M | EN-2M | ES-18UDS | ES-17UDS



ES-02JH



EN-3M



EN-2M



ES-18UDS



Specification	EN-3M	EN-2M	Features	EN-3M	EN-2M
Overall bed length	L : 2140mm	L : 2140mm	Chair Position	✓	✓
Overall bed width	W : 970mm	W : 1070mm	Foldable Type	✓	
Overall mattress base length	L : 1935mm	L : 1900mm	Assembly Type		✓
Overall mattress base width	W : 775mm	W : 910mm	Remote Controller	✓	✓
Overall bed height	4"Lockable Caster 300~700mm	4"Lockable Caster 300~700mm	Rotary Armrest Siderail	✓	✓
Head Elevation	0~70°	0~75°	3/4 Length Siderail	✓	✓
Knee Elevation	0~28°	0~15°	Plug-in Siderail	✓	✓
Trend./Rev. Trend.	11°	11°	Dust-Proof Lockable Castor	□	✓
Work Load	170kg	173.5kg	Lockable Castor	✓	
Castor Diameter	4"	4"	IV Pole Hole	✓	✓
Actuator	4	3	Mattress Stopper		✓
Volts & Frequency	AC110 ; 50/60Hz	AC110V ; 50/60Hz	Backup Battery (option) 9V Alkaline Battery x3	□	□
Actuator/Control Box	IP54	IP54			

☒ Standard
 ☐ Option


Applicable mattress size (mm)	L1980xW820xH150	L1980xW890xH100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


08 series




ES-08FDS



ES-08FSS




MS-08FS





ES-08FDS

Electric Hospital Bed



4" castor
H: 420~620
5" castor
H: 430~630




unit: mm



3/4 Length Siderail



Remote Controller




Drainage Bag Hanger



Plug-in
Head/Foot Board




Dust-Proof Central Braking
Steering Module



Mattress Stopper

ES-08FSS

Electric Hospital Bed
(Metal Frame)



Metal Frame




Remote Controller



IV Pole Hole




Lockable Castor



Bed-end Extension
(Option)

MS-08FS / MF-08FS



3 Manual Cranks / 2 Manual Cranks



4" castor
H: 420~620

5" castor
H: 430~630


MF-08FS H: 450



unit: mm




Manual Crank



Drainage Bag Hanger



3/4 Length Siderail



Mattress Stopper



Specification	ES-08FDS	ES-08FSS	MS-08FS	MF-08FS
Overall bed length	L : 2010mm	L : 1960mm	L : 2010mm	L : 2010mm
Overall bed width	W : 940mm	W : 930mm	W : 940mm	W : 940mm
Overall mattress base length	L : 1860mm	L : 1770mm	L : 1860mm	L : 1860mm
Overall mattress base width	W : 845mm	W : 815mm	W : 845mm	W : 845mm
Overall bed height	H:420~620mm (4" castor) H:430~630mm (5" castor)	H:420~620mm (4" castor) H:430~630mm (5" castor)	H:420~620mm (4" castor) H:430~630mm (5" castor)	H : 450mm
Head Elevation	0~80°	0~85°	0~80°	0~80°
Knee Elevation	0~40°	0~40°	0~40°	0~40°
Work Load	170kg	170kg	170kg	170kg
Castor Diameter	5" or 4"	5" or 4"	5" or 4"	5"
Actuator	3	3	Manual Crank	Manual Crank
Volts	AC 115V/220V	AC 115V/220V	none	none
Frequency	50/60Hz	50/60Hz	none	none
Backup Battery	none	none	none	none
Actuator/Control Box	IP42	IP42	none	none
Safety Standard	none	none	none	none

Features	ES-08FDS	ES-08FSS	MS-08FS	MF-08FS
Remote Controller	☒	☒		
3/4 Length Siderail	☒	☒	☒	☒
Dust-Proof Central Braking/Steering Module	☐	☐		
Dust-Proof Lockable Castor	☐	☐	☐	☐
Lockable Castor	☒	☒	☒	☒
Bed-Head Extension		☒		
IV Pole Hole	☒	☒	☒	☒
Drainage Bag Hanger	☒	☒	☒	☒
Mattress Stopper	☒		☒	☒

☒ Standard ☐ Option


Applicable mattress size (mm)	L1850xW840xH100	L1800xW810xH100	L1850xW840xH100	L1850xW840xH100
-------------------------------	-----------------	-----------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Mattress

JS-301, JS-302, JS-303

Mattress cover: Nylon



■ BS



Ventilation



Waterproof




Flame
Retardant



Anti
Bacterial

■ CA



Ventilation




Waterproof




Flame
Retardant


High-Density Foam



JS-301



JS-302



JS-303

Specification	JS-301		JS-302		JS-303	
Model	JS-301BS	JS-301CA	JS-302BS	JS-302CA	JS-303BS	JS-303CA
Mattress cover	■	■	■	■	■	■
Fire safety standard	BS-7177	CA-117	BS-7177	CA-117	BS-7177	CA-117
Anti-bacterial	✓		✓		✓	

JT-100

Orthopedic Traction Frame



Triangle Ring Set



Pulley Mounting Set



Foot Support Unit



Support Rod Set



1. The traction frame is suitable for all series of Joson-Care's bed and installed in four IV Pole holder.
2. The bottom of the frame is made by solid stainless steel bars.
3. Switching the wrench could quickly adjust the frame, it's tools free to set the frame.
4. The traction frame includes one set of triangular ring and pulley.
5. The octagonal pipes are made of aluminum.
6. Two sets of reinforcing units are helping to reinforce the traction frame.
7. Adjustable height and length of triangle ring set can avoid patient being impacted while position changed.

IV Stand



JR-016 IV Stand

Size : H1630~2300mm
Frame : Stainless 304#



JR-018 Stackable IV Stand

Size : 1335~2169mm
Frame : Chrome Plated steel



K5000-MS 2 Claws IV Pole

Size : 855~1510mm
Frame : Chrome Plated steel



K4000-MS 4 Claws IV Pole

Size : 1020~1760mm
Frame : Chrome Plated steel


JA-107

Emergency Cart

Overall Size : L735 x W560 x H990mm

Table Surface Size : 380 x 817 mm

Drawers size : 60 x 9 x 15 : ±1 cm




Laptop Drawer
(390 x 380mm)



Trash Bin(9 liters)
(W277 x H182 x D390mm)



Cardiac Board
(L700 x W400 x D15mm)



Disposable Lock

- Frame: Steel
- Cover: ABS
- A buckle fixed design which prevent drawers sliding out.
- Key and disposable lock available.
- Vertical and horizontal divider included.
- Flexible confuguration design.
- 2 sets divider quipe in both 3" and 6" drawers.

Accessories Option



Storage Box
(W600 x H304 x D124mm)



IV Pole
(1020~1760mm)



Suction Tube Holder
(W90 x H480 x D90mm)



Oxygen Tank Holder
(W125 x H370 x D175mm)




Pedal Trashbin(30 liters)
(W311 x H406 x D500mm)




Hand Sanitizer Holder
(Including hand sanitizer bottle)
(W75 x H105 x D77mm)



Syringe Collection Box Holder
(W130 x H150 x D90mm)



Plastic Bin Holder
(W480 x H120 x D120mm)



Smesh Bin Holder
(W488 x H125 x D125mm)



Glove Dispenser

UD Medical Cart



Medical Cart



JO-103 Mayo Table

Size : L485 x W325 x H800~1200mm
Frame : Stainless 304#



JA-108 Dressing Trolley

Size : L600 x W450 x H850mm
Frame : Stainless 304#

JA-111 Utility Cart

Size : L600 x W450 x H850mm
Frame : Stainless 304#



JF-014

Bedside Cabinet



Size : L460 x W460 x H860mm
 Color : Blue / Pink
 Material: ABS

Outpatient Equipment



JR-017 Screen

Unfolded Dimensions :
 L2100 x W470 x H1800 mm

Left and right panels :
 L600 x H1800 mm
 Center panel: L900 x H1800 mm

Materials :
 Main body made of stainless steel 304#

JO-104 Pedal

Size : Single-Layer (High) – L450 x W300 x H200mm
 Single-Layer (Low) – L450 x W300 x H150mm
 Double-Layer – L400 x W200 x H500mm
 (First step height: 230mm; Second step height: 500mm)
 Materials : Frame and legs made of stainless steel 304#



Overbed Table

JR-032 Overbed Table

Size : L815 x W403 x H820~1225mm
 Frame : Powder coated steel
 Table Surface : flame-retardant laminate
 Max. Loading : 20KG



Gas spring driven design



2" lockable castor*4
 2 with brake



JR-011 Overbed Table

Size : L765 x W385 x H720~1120mm
 Frame : Powder coated steel
 Table Surface : flame-retardant laminate

JR-051 Overbed Table

Size : L845 x W400 x H780~1050mm
 Package : L880 x W465 x H165mm
 Frame : Powder coated steel
 Table Surface : flame-retardant laminate
 Max. Loading : 20KG



Y0307-SY Overbed Table

Size : L890~1060 x W320 x H30mm
 Frame : ABS plastic
 Tube : Aluminum
 Hooks on the backside.

Doctor Chair



JO-107 Doctor Chair

Size : Ø380mm · H510~660mm
Seat: PU Foam
Frame: Aluminum



JO-109 Doctor Chair

Size : Ø300mm · H500~670mm
Frame : Stainless 304#



JO-108 Doctor Chair

Size : Ø380mm · H510~660mm
Seat: PU Foam
Frame: Aluminum

Caregiver Recliner



JD-010 Caregiver Recliner

Size :
Folded L800 x W600 x H950 mm
Unfolded L1850 x W600 x H580 mm

Materials :

- Entire structure made of stainless steel 304#
- Mattress topped with high-density foam and covered with breathable synthetic leather



CONTENTS

[illegible]



Joson-Care

We serve all your care



Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031
Release Date: 2024.06.01
Version: L